

Product Datasheet

Mouse CXCL10/IP-10 (Interferon Gamma Induced Protein 10kDa) ELISA Kit (orb3131217)

Catalog Number	orb3131217
Category	Assays and Kits
Description	Mouse CXCL10/IP-10 (Interferon Gamma Induced Protein 10kDa) ELISA Kit
Target	CXCL10
Reactivity	Mouse
UniProt ID	P17515
Assay Type	Sandwich
Assay Time	120 minutes
Range	31.25-2000pg/ml
Sensitivity	18.75pg/ml
Sample Types	Serum, Plasma, Cell Culture Supernatant, cell or tissue lysate, Other liquid samples
Storage	2 - 8 °C for 12 months
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)