

Product Datasheet

Rat GLP2 (Glucagon Like Peptide 2) ELISA Kit (orb3131238)

Catalog Number	orb3131238
Category	Assays and Kits
Description	Rat GLP2 (Glucagon Like Peptide 2) ELISA Kit
Target	GLP2
Reactivity	Rat
UniProt ID	P06883
Assay Type	Sandwich
Assay Time	90 minutes
Range	31.25-2000pg/ml
Sensitivity	18.75pg/ml
Sample Types	Serum, Plasma, Cell Culture Supernatant, cell or tissue lysate, Other liquid samples
Storage	2 - 8 °C for 12 months
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)