

Product Datasheet

Others,PD-1/PD-L1 Antibody (orb3131351)

Catalog Number	orb3131351
Category	Antibodies
Description	IBI-322 is a bispecific antibody targeting PD-L1 and CD47. In monovalent binding, it can attenuate the activity of CD47, and in bivalent binding, it can block the activity of PD-L1. IBI-322 selectively binds to CD47+PD-L1+ tumor cells, effectively inhibiting the CD47-SIRPα signaling and inducing strong tumor cell phagocytosis in vitro.
Target	Others,PD-1/PD-L1
Clonality	Recombinant
Conjugation	Unconjugated
Antibody Type	Inhibitor Antibody
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558