

Product Datasheet

TGF-beta/Smad Antibody (orb3131495)

Catalog Number	orb3131495
Category	Antibodies
Description	SAR-439459 is a human monoclonal antibody (mAb) targeting TGFβ1/TGFβ-1. It blocks TGFβ-mediated pSMAD signaling and inhibits T cell and NK cell activity. Additionally, SAR-439459 exhibits antitumor properties.
Target	TGF-beta/Smad
Clonality	Monoclonal
Conjugation	Unconjugated
Antibody Type	Inhibitor Antibody
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)