

Product Datasheet

AC Antibody purification resin 2 (orb3137167)

Catalog Number	orb3137167
Category	Small Molecules
Description	AC Antibody Purification Resin 2 is composed of spherical, highly cross-linked agarose beads characterized by a narrow size distribution and strong mechanical stability. It is utilized for the separation and purification of complex antibodies such as monoclonal, bispecific, and polyclonal antibodies as well as Fc fusion proteins. [Substrate: Highly cross-linked agarose microspheres; Particle size: 65µm; Ligand: Alkali-resistant recombinant Protein A; ADC Purification Resin.]
Target	Others
Storage	-20°C
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)