

Product Datasheet

di-DTPA-LTL (orb3143391)

Catalog Number	orb3143391
Category	Small Molecules
Description	di-DTPA-LTL is a divalent hapten designed based on tyrosine-containing polypeptides. It exhibits excellent hydrophilicity and biodistribution properties. di-DTPA-LTL is labeled with ¹¹¹ In (indium) and ¹³¹ I (iodine). By administering a fixed low dose (5 mg) of a bispecific antibody (anti-CEA x, anti-DTPA) and the di-DTPA hapten peptide (labeled with ¹¹¹ In) to patients, di-DTPA-LTL enables tumor radioimmunoimaging in patients with primary colorectal cancer expressing CEA.
Target	Others
Formula	C ₅₂ H ₈₂ N ₁₂ O ₂₃
Storage	-20°C
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)