

Product Datasheet

Mouse VEGF121 Protein, His tag (Animal-Free) (orb3148457)

Catalog Number	orb3148457
Category	Proteins
Description	Vascular Endothelial Growth Factors 121 (VEGF121) is a truncated version of VEGF165, which produced in E. coli is a homodimer, non-glycosylated, polypeptide chain and having a molecular mass of 28.4 kDa. There is three different isoforms (120, 164 and 188 a.a.) found in mouse. VEGF 121 shows that lack basic heparin-binding regions and are freely diffusible. Mouse VEGF121 shares 98% identity with corresponding regions of rat, 89% with canine, feline, equine and porcine, and 87% with human, ovine and bovine VEGF, respectively.
Storage	Please store the whole kit at 4°C upon receipt ,and protect from light. The kit is valid for 28 months. Please refer to the table below Materials supplied and Storage conditions to store all the components.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)