

Product Datasheet

Human TNF- β Protein, His tag (Animal-Free) (orb3148509)

Catalog Number	orb3148509
Category	Proteins
Description	Human TNF beta (Lymphotoxin-alpha) is a 22 kDa cytokine with 205 amino acid residues. TNF beta belongs to tumor necrosis factor (TNF) superfamily as TNF alpha. Binding to lymphotoxin- β receptor (Lt β R), canonical and non-canonical NF- κ B signal pathway will be activated. TNF beta is related to apoptosis and inflammatory signals
Storage	Lyophilized protein product should be stored desiccated below -18°C. Upon reconstitution, the protein should be stored at 4°C between 2-7 days and for future use below -18°C. For long term storage it is recommended to add a carrier protein (5% HSA,10% FBS or 0.1%BSA). Please prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)