

Product Datasheet

Human MDA5/IFIH1 Antibody (orb3149164)

Catalog Number	orb3149164
Category	Antibodies
Description	Anti-Human MDA5/IFIH1 Antibody (SAA2314) [SAA2314] is a human monoclonal antibody detecting Interferon-induced helicase C domain-containing protein 1 in ELISA, WB . Suitable for Human .
Species/Host	Human
Conjugation	Unconjugated
Reactivity	Human
Form/Appearance	Liquid
Purity	>95% as determined by SDS-PAGE.
UniProt ID	Q9BYX4
Tested applications	ELISA, WB
Antibody Type	Primary Antibody
Clone Number	SAA2314
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)