

Product Datasheet

Recombinant Mouse CD40/TNFRSF5 Protein, C-Fc (orb3149219)

Catalog Number	orb3149219
Category	Proteins
Description	Recombinant Mouse CD40/TNFRSF5 Protein, C-Fc is produced in Mammalian Cells and supplied in lyophilized form. The protein includes the regions of Met1-Arg193 and carries the C-Fc tag. With a molecular weight of 47.83 kDa, it is suitable for ELISA, SDS-PAGE and WB applications. This product has a purity of greater than 95% as determined by SDS-PAGE.
Form/Appearance	Lyophilized
Purity	>95% as determined by SDS-PAGE.
UniProt ID	P27512
MW	47.83 kDa
Tested applications	ELISA, SDS-PAGE, WB
Expression System	Mammalian Cells
Expression Region	Met1-Arg193
Storage	-20°C
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.


7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558



SDS-PAGE for Recombinant Mouse CD40/TNFRSF5 Protein, C-Fc.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)