

## Product Datasheet

# Recombinant Human UBE2H Protein, N-His (orb3150064)

<b>Catalog Number</b>	orb3150064
<b>Category</b>	Proteins
<b>Description</b>	Recombinant Human UBE2H Protein, N-His expressed in E. coli, spanning Met1-Leu183. Purity: >90% by SDS-PAGE. Highlights His-Tagged: N-terminal 6×His tag for IMAC purification. E. coli Expression: High-yield, cost-effective production. High Purity: >90% purity verified by SDS-PAGE.
<b>Form/Appearance</b>	Lyophilized
<b>Purity</b>	>90% as determined by SDS-PAGE.
<b>UniProt ID</b>	<b>P62256</b>
<b>MW</b>	22.82 kDa
<b>Tested applications</b>	ELISA, SDS-PAGE, WB
<b>Expression System</b>	E. coli
<b>Expression Region</b>	Met1-Leu183
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Expiration Date</b>	6 months from date of receipt.

### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)