

## Product Datasheet

### FIP1L1 Polyclonal Antibody (orb3150399)

|                            |  |
|----------------------------|--|
| <b>Catalog Number</b>      | orb3150399   |
| <b>Category</b>            | Antibodies   |
| <b>Description</b>         | Anti-FIP1L1 Polyclonal Antibody is a rabbit polyclonal antibody detecting Pre-mRNA 3'-end-processing factor FIP1; hFip1; FIP1-like... in ELISA, IHC, WB. Suitable for Human, Rat, and Mouse. Highlights Affinity Purified: Minimal background and high purity for reliable results. Multi-Application: Validated across multiple applications. Multi-Species: Cross-reactive for translational research. |
| <b>Clonality</b>           | Polyclonal   |
| <b>Species/Host</b>        | Rabbit   |
| <b>Conjugation</b>         | Unconjugated   |
| <b>Reactivity</b>          | Human, Mouse, Rat  |
| <b>Form/Appearance</b>     | Liquid   |
| <b>Concentration</b>       | Lot dependent  |
| <b>UniProt ID</b>          | <b>Q5U317, Q9D824, Q6UN15</b>  |
| <b>Tested applications</b> | ELISA, IHC, WB   |
| <b>Antibody Type</b>       | Primary Antibody   |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.  |
| <b>Note</b>                | For research use only  |
| <b>Expiration Date</b>     | 12 months from date of receipt.  |

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)