

Product Datasheet

CHIKV Spike glycoprotein E2 Antibody (orb3151115)

Catalog Number	orb3151115
Category	Antibodies
Description	Anti-CHIKV Spike glycoprotein E2 Antibody (BE59) is a human monoclonal antibody detecting p62 in ELISA. Suitable for Chikungunya virus.
Species/Host	Human
Conjugation	Unconjugated
Reactivity	Virus
Form/Appearance	Liquid
Concentration	Lot dependent
Purity	>95% as determined by SDS-PAGE.
UniProt ID	Q8JUX5
Tested applications	ELISA
Antibody Type	Primary Antibody
Clone Number	BE59
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.


7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558



SDS-PAGE for Anti-CHIKV Spike glycoprotein E2 Antibody.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)