

Product Datasheet

Human IL1RAP & CD3 Bispecific Antibody (orb3151179)

| | |
|----------------------------|--|
| Catalog Number | orb3151179 |
| Category | Antibodies |
| Description | Anti-Human IL1RAP & CD3 Bispecific Antibody (BDKH224) is a human monoclonal antibody detecting Interleukin-1 receptor accessory protein in ELISA, FCM, in vivo study. Suitable for Human. Highlights Multi-Application: Validated across multiple applications. Flow Cytometry: Validated for flow cytometry applications. |
| Clonality | Recombinant |
| Species/Host | Human |
| Conjugation | Unconjugated |
| Reactivity | Human |
| Form/Appearance | Liquid |
| Concentration | Lot dependent |
| Purity | >95% as determined by SDS-PAGE. |
| UniProt ID | Q9NPH3, P07766 |
| Tested applications | ELISA, FC |
| Antibody Type | Recombinant Antibody |
| Clone Number | BDKH224 |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Storage


Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note

For research use only

Expiration Date

12 months from date of receipt.



SDS-PAGE for Anti-Human IL1RAP & CD3 Bispecific Antibody.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)