

Product Datasheet

Mouse PD1 / VEGF-A Bispecific Antibody (orb3152001)

Catalog Number	orb3152001
Category	Antibodies
Description	Anti-Mouse PD1 / VEGF-A Bispecific Antibody (B41HX3C34XI8X53) is a mouse monoclonal antibody detecting Programmed cell death protein 1 in ELISA, FACS, Functional assay, Research in vivo. Suitable for Mouse. Highlights Multi-Application: Validated across multiple applications.
Clonality	Recombinant
Species/Host	Mouse
Conjugation	Unconjugated
Reactivity	Mouse
Form/Appearance	Liquid
Concentration	Lot dependent
Purity	>95% as determined by SDS-PAGE.
UniProt ID	Q00731 , Q02242
Tested applications	ELISA, FA, FACS, In vivo
Antibody Type	Primary Antibody
Clone Number	B41HX3C34XI8X53
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC


68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Note For research use only

Expiration Date 12 months from date of receipt.



SDS-PAGE for Anti-Mouse PD1 / VEGF-A Bispecific Antibody.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)