

Product Datasheet

PBP2a Antibody (orb3152051)

Catalog Number	orb3152051
Category	Antibodies
Description	Anti-PBP2a Antibody (BUCC4639) is a rabbit monoclonal antibody detecting Beta-lactam-inducible penicillin-binding protein in ELISA, WB, SPR. Highlights Rabbit Monoclonal: High specificity, affinity, and excellent lot-to-lot consistency. Multi-Application: Validated across multiple applications.
Species/Host	Rabbit
Conjugation	Unconjugated
Reactivity	Bacteria
Form/Appearance	Liquid
Concentration	Lot dependent
Purity	>95% as determined by SDS-PAGE.
UniProt ID	P07944
Tested applications	ELISA, SBR, WB
Antibody Type	Primary Antibody
Clone Number	BUCC4639
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com


Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)



SDS-PAGE for Anti-PBP2a Antibody.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)