

## Product Datasheet

### PBP2a Antibody (orb3152052)

<b>Catalog Number</b>	orb3152052
<b>Category</b>	Antibodies
<b>Description</b>	Anti-PBP2a Antibody (SAA2418) [SAA2418] is a rabbit monoclonal antibody detecting Beta-lactam-inducible penicillin-binding protein in ELISA, WB, SPR. Highlights Rabbit Monoclonal: High specificity, affinity, and excellent lot-to-lot consistency. Multi-Application: Validated across multiple applications.
<b>Species/Host</b>	Rabbit
<b>Conjugation</b>	Unconjugated
<b>Reactivity</b>	Bacteria
<b>Form/Appearance</b>	Liquid
<b>Purity</b>	>95% as determined by SDS-PAGE.
<b>UniProt ID</b>	<b>P07944</b>
<b>Tested applications</b>	ELISA, SBR, WB
<b>Antibody Type</b>	Primary Antibody
<b>Clone Number</b>	SAA2418
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Expiration Date</b>	12 months from date of receipt.

#### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

#### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)