

Product Datasheet

Human CD41/ITGA2B & CD61/ITGB3 Antibody, PerCP (orb3152089)

Catalog Number	orb3152089
Category	Antibodies
Description	Anti-Human CD41/ITGA2B & CD61/ITGB3 Antibody (PAC-1), PerCP [PAC-1] is a mouse monoclonal antibody detecting CD41 in FCM. Suitable for Human. Highlights Flow Cytometry — Validated for flow cytometry applications. Fluorescent Label — Direct conjugate for streamlined workflows.
Species/Host	Mouse
Conjugation	PerCP
Reactivity	Human
Form/Appearance	Liquid
Purity	>95% as determined by SDS-PAGE.
UniProt ID	P05106, P08514
Tested applications	FC
Clone Number	PAC-1
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)