

Product Datasheet

Human CD2AP Monoclonal Antibody (orb3152725)

| | |
|----------------------------|---|
| Catalog Number | orb3152725 |
| Category | Antibodies |
| Description | Anti-Human CD2AP Monoclonal Antibody (B3C897) is a mouse monoclonal antibody detecting CD2-associated protein in ELISA, IHC, WB. Suitable for Human. Highlights Recombinant: Superior batch-to-batch reproducibility with defined sequences. Multi-Application: Validated across multiple applications. |
| Clonality | Monoclonal |
| Species/Host | Mouse |
| Conjugation | Unconjugated |
| Reactivity | Human |
| Form/Appearance | Liquid |
| UniProt ID | Q9Y5K6 |
| Tested applications | ELISA, IHC, WB |
| Antibody Type | Primary Antibody |
| Clone Number | B3C897 |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Expiration Date | 12 months from date of receipt. |

Biorbyt Ltd.


7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)



SDS-PAGE for Anti-Human CD2AP Monoclonal Antibody.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)