

Product Datasheet

Recombinant Human CD254/RANKL/TNFSF11 Protein, N-His (orb3154534)

Catalog Number	orb3154534
Category	Proteins
Description	Recombinant Human CD254/RANKL/TNFSF11 Protein, N-His expressed in Mammalian Cells, spanning Gly136-Asp317. Purity: >90% by SDS-PAGE. Highlights His-Tagged: N-terminal 6×His tag for IMAC purification. Mammalian Expression: Proper folding and native glycosylation. High Purity: >90% purity verified by SDS-PAGE.
Form/Appearance	Lyophilized
Purity	>90% as determined by SDS-PAGE.
UniProt ID	O14788
MW	26.92 kDa
Tested applications	ELISA, SDS-PAGE, WB
Expression System	Mammalian Cells
Expression Region	Gly136-Asp317
Storage	-20°C
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)