

Product Datasheet

Recombinant Mouse KLK7 Protein, C-His (orb3154836)

Catalog Number	orb3154836
Category	Proteins
Description	Recombinant Mouse KLK7 Protein, C-His is produced in Mammalian Cells and supplied in lyophilized form. The protein includes the regions of Gln22-Arg249 and carries the C-His tag. With a molecular weight of 27.98 kDa, it is suitable for ELISA, SDS-PAGE and WB applications. This product has a purity of greater than 90% as determined by SDS-PAGE.
Form/Appearance	Lyophilized
Purity	>90% as determined by SDS-PAGE.
UniProt ID	Q91VE3
MW	27.98 kDa
Tested applications	ELISA, SDS-PAGE, WB
Expression System	Mammalian Cells
Expression Region	Gln22-Arg249
Storage	-20°C
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)