

## Product Datasheet

### Recombinant Mouse KLK7 Protein, C-Fc (orb3154837)

|                            |                                                                                                                                                                                                                                                                                                                                                         |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog Number</b>      | orb3154837                                                                                                                                                                                                                                                                                                                                              |
| <b>Category</b>            | Proteins                                                                                                                                                                                                                                                                                                                                                |
| <b>Description</b>         | Recombinant Mouse KLK7 Protein, C-Fc is produced in Mammalian Cells and supplied in lyophilized form. The protein includes the regions of Gln22-Arg249 and carries the C-Fc tag. With a molecular weight of 53.41 kDa, it is suitable for ELISA, SDS-PAGE and WB applications. This product has a purity of greater than 90% as determined by SDS-PAGE. |
| <b>Form/Appearance</b>     | Lyophilized                                                                                                                                                                                                                                                                                                                                             |
| <b>Purity</b>              | >90% as determined by SDS-PAGE.                                                                                                                                                                                                                                                                                                                         |
| <b>UniProt ID</b>          | <b>Q91VE3</b>                                                                                                                                                                                                                                                                                                                                           |
| <b>MW</b>                  | 53.41 kDa                                                                                                                                                                                                                                                                                                                                               |
| <b>Tested applications</b> | ELISA, SDS-PAGE, WB                                                                                                                                                                                                                                                                                                                                     |
| <b>Expression System</b>   | Mammalian Cells                                                                                                                                                                                                                                                                                                                                         |
| <b>Expression Region</b>   | Gln22-Arg249                                                                                                                                                                                                                                                                                                                                            |
| <b>Storage</b>             | -20°C                                                                                                                                                                                                                                                                                                                                                   |
| <b>Note</b>                | For research use only                                                                                                                                                                                                                                                                                                                                   |
| <b>Expiration Date</b>     | 6 months from date of receipt.                                                                                                                                                                                                                                                                                                                          |

#### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

#### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558