

## Product Datasheet

### IgG control Antibody (APC) (orb3158135)

<b>Catalog Number</b>	orb3158135
<b>Category</b>	Antibodies
<b>Description</b>	Monoclonal Rabbit IgG Antibody Isotype Control is a mouse monoclonal antibody recombinantly expressed from human 293 cells (HEK293), which provides higher batch consistency and long term security of supply.
<b>Target</b>	IgG control
<b>Species/Host</b>	Human
<b>Conjugation</b>	APC
<b>Tag</b>	Rabbit IgG   Rabbit Kappa
<b>Reactivity</b>	Rabbit
<b>Form/Appearance</b>	Powder
<b>Buffer/Preservatives</b>	PBS, pH 7.4, 0.03% Proclin300, 0.5% rAlbumin
<b>Storage</b>	-20°C, avoid light
<b>Note</b>	For research use only
<b>Expiration Date</b>	12 months from date of receipt.

#### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

#### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)