

Product Datasheet

Goat Anti-human IgG (H+L) Antibody (HRP) (orb3158275)

Catalog Number	orb3158275
Category	Antibodies
Description	Goat Anti-human IgG (H+L) Antibody (HRP) is an HRP conjugated polyclonal antibody. Generated against a Human IgG (H+L), it is supplied in liquid form and exhibits reactivity with Human samples. The suitable applications for this antibody are ELISA, IHC-P, IP and WB.
Clonality	Polyclonal
Species/Host	Goat
Isotype	IgG
Conjugation	HRP
Reactivity	Human
Form/Appearance	Liquid
Concentration	0.6 mg/ml
Buffer/Preservatives	PBS/ glycerol
Purification	Affinity-purified
Immunogen	Human IgG (H+L)
Tested applications	ELISA, IHC-P, IP, WB
Dilution range	WB:1:500-20000, IHC-P:1:200-2000, IP:1:100-1:1000, ELISA: 1:2000-50000

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Antibody Type	Secondary Antibody
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)