

## Product Datasheet

### FOXD3 Antibody (orb3160681)

|                             |   |
|-----------------------------|---|
| <b>Catalog Number</b>       | orb3160681  |
| <b>Category</b>             | Antibodies  |
| <b>Description</b>          | FOXD3 Antibody                                    |
| <b>Target</b>               | FOXD3   |
| <b>Clonality</b>            | Polyclonal  |
| <b>Species/Host</b>         | Rabbit  |
| <b>Isotype</b>              | IgG   |
| <b>Conjugation</b>          | Unconjugated                                      |
| <b>Reactivity</b>           | Human, Mouse, Rat                                 |
| <b>Form/Appearance</b>      | Liquid  |
| <b>Concentration</b>        | 1mg/ml  |
| <b>Buffer/Preservatives</b> | Purified antibody in PBS with 0.05% sodium azide. |
| <b>Purification</b>         | Protein A   |
| <b>Immunogen</b>            | Polypeptide                                       |
| <b>UniProt ID</b>           | <b>Q9UJU5</b>                                     |
| <b>MW</b>                   | 48KD  |
| <b>Tested applications</b>  | ELISA, IHC, WB                                    |
| <b>Dilution range</b>       | WB;:1:1000-1:2000;IHC;:1:200-1:2000               |

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

|                        |  |
|------------------------|--|
| <b>Antibody Type</b>   | Primary Antibody   |
| <b>DNA Notes</b>       | This gene belongs to the forkhead family of transcription factors which is characterized by a distinct forkhead domain. Mutations in this gene cause autoimmune susceptibility 1. [provided by RefSeq, Nov 2008] |
| <b>Storage</b>         | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.  |
| <b>Note</b>            | For research use only  |
| <b>Expiration Date</b> | 12 months from date of receipt.  |

---

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)