

Product Datasheet

Recombinant human IL-17A protein (Active, CHO), BF647 conjugated (orb3164953)

Catalog Number	orb3164953
Category	Proteins
Description	Recombinant human IL-17A protein (Active, CHO), BF647 conjugated
Target	Interleukin-17A
Tag	Tag free
Form/Appearance	Lyophilized or Liquid
Buffer/Preservatives	Lyophilized from a 0.22µm-filtered solution containing PBS(pH7.4) with 5% Trehalose.
Reconstitution	It is recommended to redissolve in sterile deionized water
Purity	≥ 95% as determined by SDS-PAGE.
Purification	AC
Protein Sequence	20-155/155
MW	15.5 kDa
Expression System	CHO cell
Biological Origin	Human
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. 36 months at -20°C to -80°C in lyophilized state. 6 months at -20°C to -80°C under sterile conditions after reconstitution. 7-10 days at 2°C to 35°C under sterile conditions after reconstitution.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)