

## Product Datasheet

### CRP Monoclonal Antibody (Detection) (orb316745)

|                             |   |
|-----------------------------|---|
| <b>Catalog Number</b>       | orb316745   |
| <b>Category</b>             | Antibodies  |
| <b>Description</b>          | CRP Monoclonal Antibody (Detection)   |
| <b>Target</b>               | C-Reactive Protein (CRP)  |
| <b>Clonality</b>            | Monoclonal  |
| <b>Conjugation</b>          | Unconjugated  |
| <b>Buffer/Preservatives</b> | NaN3  |
| <b>Tested applications</b>  | EIA   |
| <b>Matched Pairs</b>        | Capture: orb316746; Detection: orb316745;   |
| <b>Storage</b>              | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>                 | For research use only   |
| <b>Expiration Date</b>      | 12 months from date of receipt.   |

#### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

#### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)