

Product Datasheet

Endomorphin-2 (E-2) Antibody (orb3167625)

Catalog Number	orb3167625
Category	Antibodies
Description	Our Anti-Endomorphin-2 (E-2) rabbit polyclonal primary antibody detects human, mouse, and rat Endomorphin-2 (E-2), and is affinity-purified. It is validated for use in ICC, IHC-Frozen, IHC-Paraffin-embedded.
Target	Endomorphin-2 (E-2)
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	Unconjugated
Reactivity	Human, Mouse, Rat
Form/Appearance	Lyophilized
Buffer/Preservatives	Lyophilized with rAlbumin.
Immunogen	Synthetic Endomorphin-2 (YPFF) peptide conjugated to BSA. Selectivity towards Endomorphin-2 was achieved by separating fraction of affinity purified antibodies which did not bind to Endomorphin-1.
Tested applications	ICC, IHC
Dilution range	IHC: 5-15 µg/mL; ICC: 5-15 µg/mL
Specificity	Endomorphin-2 peptide.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)