

Product Datasheet

Annexin V-FITC Apoptosis Detection Kit (orb322264)

Catalog Number	orb322264
Category	Assays and Kits
Description	This Annexin V-FITC Apoptosis Detection Kit enables detection of apoptosis and membrane permeability. It uses the green fluorescent reagent Annexin V-FITC to label exposed phosphatidylserine on the surface of cultured cells. Cells with compromised membranes-a hallmark of late apoptosis or necrosis-are identified using the live/dead stain propidium iodide.
Assay Type	Detection
Storage	Storage Temp: 4°C
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)