

Product Datasheet

NOP5 Rabbit Polyclonal Antibody (orb326124)

Catalog Number	orb326124
Category	Antibodies
Description	NOP5 Rabbit Polyclonal Antibody is an unconjugated, affinity purified antibody for the detection of NOP58. It is suitable for WB application. The immunogen is a synthetic peptide directed towards the middle region of human NOP5/NOP58. This antibody reacts with Human samples and is predicted to react with Bovine, Canine, Equine, Mouse, Porcine, Rabbit, Rat. It is supplied as a liquid in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Target	NOP58
Clonality	Polyclonal
Species/Host	Rabbit
Conjugation	Unconjugated
Reactivity	Human
Predicted Reactivity	Bovine, Canine, Equine, Mouse, Porcine, Rabbit, Rat
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Buffer/Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Purification	Affinity Purified
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human NOP5/NOP58

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom


Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Protein Sequence	Synthetic peptide located within the following region: LRTLEDRGIRKISGTGKALAKTEKYEHKSEVKTYDPSGDSTLPTCSKKRK
UniProt ID	Q9Y2X3
MW	59kDa
Tested applications	WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
NCBI	NP_057018
Expiration Date	12 months from date of receipt.



WB Suggested Anti-NOP5 Antibody Titration: 0.2-1 ug/mL, ELISA
Titer: 1:312500, Positive Control: MCF7 cell lysate.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)