

Product Datasheet

PDK1 Rabbit Polyclonal Antibody (orb329620)

Description	Rabbit polyclonal antibody to PDK1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human PDK1
Target	PDK1
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: ATMEHHANRGVYPPIQVHVTLGNEDLTVKMSDRGGGVPLRKIDRLFNYMY
Clonality	Polyclonal
MW	47kDa
Uniprot ID	Q15118

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.


68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

NCBI

NP_002601


Expiration Date

12 months from date of receipt.




Host: Rabbit
Target Name: Pdk1
Sample Type: Mouse Testis Lysate
Antibody Dilution: 1.0µg/ml

Sample Tissue: Mouse Testis, Antibody Dilution: 1 ug/mL.



Rabbit Anti-PDK2
Sample Type: Human Adult Placenta
Antibody Concentration: 1ug/mL

Sample Type: Human Adult Placenta, Antibody Dilution: 1.0 ug/mL.



Host: Rabbit
Target Name: Pdk1
Sample Type: Mouse Testis Lysate
Antibody Dilution: 1.0µg/ml


Sample Tissue: Mouse Testis, Antibody Dilution: 1 ug/mL.

Biorbyt Ltd.


7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)



Sample Type: Human Jurkat, Antibody Dilution: 1.0 ug/mL, PDK1 is supported by BioGPS gene expression data to be expressed in Jurkat.



WB Suggested Anti-PDK1 Antibody Titration: 0.2-1 ug/mL, ELISA Titer: 1:1562500, Positive Control: MCF7 cell lysate.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)