

## Product Datasheet

### DDX6 Rabbit Polyclonal Antibody (orb329867)

|                             |   |
|-----------------------------|---|
| <b>Catalog Number</b>       | orb329867   |
| <b>Category</b>             | Antibodies  |
| <b>Description</b>          | Rabbit polyclonal antibody to DDX6  |
| <b>Target</b>               | DDX6  |
| <b>Clonality</b>            | Polyclonal  |
| <b>Species/Host</b>         | Rabbit  |
| <b>Conjugation</b>          | Unconjugated  |
| <b>Reactivity</b>           | Human   |
| <b>Predicted Reactivity</b> | Bovine, Canine, Equine, Goat, Guinea pig, Mouse, Rabbit, Rat, Yeast, Zebrafish                              |
| <b>Form/Appearance</b>      | Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.           |
| <b>Concentration</b>        | 0.5 mg/ml   |
| <b>Buffer/Preservatives</b> | Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.           |
| <b>Purification</b>         | Affinity Purified   |
| <b>Immunogen</b>            | The immunogen is a synthetic peptide directed towards the N terminal region of human DDX6                   |
| <b>Protein Sequence</b>     | Synthetic peptide located within the following region:<br>CLKRELLMGIFEMGWKPSPIQEESIPIALSGRDILARAKNGTGKSGAYL |
| <b>UniProt ID</b>           | <b>P26196</b>   |

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)


**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

|                            |   |
|----------------------------|---|
| <b>MW</b>                  | 54kDa   |
| <b>Tested applications</b> | WB  |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>                | For research use only   |
| <b>NCBI</b>                | <b>NP_004388</b>  |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |

**DDX6**

Sample Tissue: Human Hela, Antibody Dilution: 1.0 ug/mL.



WB Suggested Anti-DDX6 Antibody Titration: 0.2-1 ug/mL, ELISA Titer: 1:312500, Positive Control: OVCAR-3 cell lysate, DDX6 is supported by BioGPS gene expression data to be expressed in OVCAR3.

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)