

Product Datasheet

Ppargc1a Rabbit Polyclonal Antibody (orb330035)

Catalog Number	orb330035
Category	Antibodies
Description	Rabbit polyclonal antibody to Ppargc1a
Target	Ppargc1a
Clonality	Polyclonal
Species/Host	Rabbit
Conjugation	Unconjugated
Reactivity	Human, Mouse
Predicted Reactivity	Bovine, Canine, Equine, Goat, Guinea pig, Rabbit, Rat, Zebrafish
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Buffer/Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Purification	Affinity Purified
Immunogen	The immunogen is a synthetic peptide corresponding to a region of Mouse
Protein Sequence	Synthetic peptide located within the following region: LENGYTLRRSNETDFELYFCGRKQFFKSNYADLDTNSDDFDPASTKSKYD
UniProt ID	O70343

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)


MW 90kDa

Tested applications WB

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

NCBI [NP_032930](#)



Lanes: Lane 1: 20 ug cytosolic fraction mesangial cells, Lane 2: 20 ug cytosolic fraction mesangial cells, Lane 3: 20 ug mitochondrial fraction mesangial cells, Lane 4: 20 ug treated cytosolic fraction mesangial cells, Lane 5: 20 ug Myometrial tissue lysate, Primary Antibody Dilution: 1:500, Secondary Antibody: Goat anti-rabbit HRP, Secondary Antibody Dilution: 1:10000, Gene Name: PGC1 alpha.



WB Suggested Anti-Ppargc1α Antibody Titration: 0.2-1 ug/mL, ELISA Titer: 1:312500, Positive Control: Mouse Brain.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)