

Product Datasheet

PERP Rabbit Polyclonal Antibody (orb330597)

Description	Rabbit polyclonal antibody to PERP
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human PERP
Target	PERP
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: FLRVIGLLALAAVFQIISLVIYPVKYTQTFTLHANPAVTYIYNWAYGFG
Clonality	Polyclonal
MW	21kDa
Uniprot ID	Q96FX8

Biorbyt Ltd.


7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+441223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

NCBI**NP_071404****Expiration Date**

12 months from date of receipt.



Rabbit Anti-PERP Antibody, Catalog Number: orb330597,
Formalin Fixed Paraffin Embedded Tissue: Human heart Tissue,
Observed Staining: Plasma membrane in intercalated disk,
Primary Antibody Concentration: 1:100, Other Working
Concentrations: 1:600, Secondary Antibody: Donkey anti-
Rabbit-Cy3, Secondary Antibody Concentration: 1:200,
Magnification: 20X, Exposure Time: 0.5-2.0 sec.



WB Suggested Anti-PERP Antibody Titration: 0.2-1 ug/ml, ELISA
Titer: 1:62500, Positive Control: PANC1 cell lysate, PERP is
supported by BioGPS gene expression data to be expressed in
PANC1.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)