

Product Datasheet

Anti-FOLR2 Purified (orb334037)

Catalog Number	orb334037
Description	Mouse monoclonal antibody to FOLR2
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, IP, WB
Immunogen	BW5147alpha,beta- cells
Target	FOLR2
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Western blotting: Non-reducing conditions. Flow cytometry: Recommended dilution: 1-4 µg/ml
Isotype	Mouse IgG1
Clonality	Monoclonal
Clone Number	EM-35
Antibody Type	Primary Antibody
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P14207

Biorbyt Ltd.


7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.


68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Entrez**2350****Expiration Date**

12 months from date of receipt.



Western blotting analysis of FOLR2 in THP1 cells (A) and FOLR2-transfected THP1 cells (B) using anti-FOLR2 (EM-35) purified.



Separation of FOLR2 transfected BW5147 cells stained using anti-FOLR2 (EM-35) purified antibody (concentration in sample 4 µg/ml, GAM APC, red-filled) from FOLR2 transfected BW5147 cells unstained by primary antibody (GAM APC, black-dashed) in flow cytometry analysis (surface staining).

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)