

Product Datasheet

Hepatitis B Virus Surface Antigen Goat Polyclonal Antibody (Biotin) (orb334871)

Catalog Number	orb334871
Category	Antibodies
Description	The Hepatitis B Virus Surface Antigen Antibody-Biotin labeled is suitable for ELISA. It is a Polyclonal, Biotin conjugated antibody which raised against Recombinant Hepatitis B Virus Surface Antigen. Purification: Immunogen affinity purified
Clonality	Polyclonal
Species/Host	Goat
Conjugation	Biotin
Reactivity	Virus
Form/Appearance	Liquid
Buffer/Preservatives	0.01M Phosphate Buffered Saline, pH 7.2, 50% glycerol, 0.05% Sodium Azide.
Purification	Immunogen affinity purified
Immunogen	Recombinant Hepatitis B Virus Surface Antigen
Tested applications	ELISA
Dilution range	ELISA: Use at an assay dependent dilution.
Specificity	Reacts only with Hepatitis B Virus Surface Antigen.
Antibody Type	Primary Antibody

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)