

Product Datasheet

Human lipoprotein ELISA Kit (orb339715)

Catalog Number	orb339715
Category	Assays and Kits
Description	The Human lipoprotein ELISA Kit detects the target LPA in human samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 62.5 ng/mL-4000 ng/mL. The kit offers a sensitivity of 53.106 ng/mL.
Target	LPA
Reactivity	Human
Application notes	Detection Wavelength: 450 nm
Assay Type	Competitive
Assay Time	1-5h
Range	62.5 ng/mL-4000 ng/mL
Sensitivity	53.106 ng/mL
Sample Types	serum, plasma, tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)