

Product Datasheet

Mouse DKK1 ELISA Kit (orb339739)

Catalog Number	orb339739
Category	Assays and Kits
Description	The Mouse DKK1 ELISA Kit detects the target DKK1 in mouse samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 78.125 pg/mL-5000 pg/mL. The kit offers a sensitivity of 13.127 pg/mL.
Target	DKK1
Reactivity	Mouse
UniProt ID	O54908
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	78.125 pg/mL-5000 pg/mL
Sensitivity	13.127 pg/mL
Sample Types	serum, plasma, tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)