

## Product Datasheet

### Human MT ELISA Kit (orb339769)

<b>Catalog Number</b>	orb339769
<b>Category</b>	Assays and Kits
<b>Description</b>	The Human MT ELISA Kit detects the target MT in human samples. It measures levels in sample types including serum, plasma, cell culture supernates and tissue homogenates and urine and saliva. The detection range is 6.25 pg/mL-400 pg/mL. The kit offers a sensitivity of 0.813 pg/mL.
<b>Target</b>	MT
<b>Reactivity</b>	Human
<b>UniProt ID</b>	<b>P04731</b>
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Assay Type</b>	Sandwich
<b>Assay Time</b>	1-5h
<b>Range</b>	6.25 pg/mL-400 pg/mL
<b>Sensitivity</b>	0.813 pg/mL
<b>Sample Types</b>	serum, plasma, cell culture supernates, tissue homogenates, urine, saliva
<b>Note</b>	For research use only
<b>Expiration Date</b>	Please enquire.

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)