

Product Datasheet

Mouse CD48 ELISA Kit (orb339821)

Catalog Number	orb339821
Category	Assays and Kits
Description	The Mouse CD48 ELISA Kit detects the target CD48 in mouse samples. Assay time: 1-5h.
Target	CD48
Reactivity	Mouse
UniProt ID	P18181
Application notes	Detection Wavelength: 450 nm
Assay Time	1-5h
Range	Request Information
Sensitivity	Request Information
Sample Types	Request Information
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)