

Product Datasheet

Pig TNF/TNFA/TNFSF2 ELISA Kit (orb339833)

Catalog Number	orb339833
Category	Assays and Kits
Description	The Pig TNF/TNFA/TNFSF2 ELISA Kit detects the target TNF in porcine samples. It measures levels in sample types including serum, plasma, cell culture supernates and tissue homogenates. The detection range is 0.047 ng/mL-3 ng/mL. The kit offers a sensitivity of 0.040 ng/mL.
Target	TNF
Reactivity	Porcine
UniProt ID	P23563
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	0.047 ng/mL-3 ng/mL
Sensitivity	0.040 ng/mL
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)