

Product Datasheet

Bovine LGB ELISA Kit (orb339837)

Catalog Number	orb339837
Category	Assays and Kits
Description	The Bovine LGB ELISA Kit detects the target LGB in bovine samples. It measures levels in sample types including serum and plasma and milk. The detection range is 0.156 µg/mL-10 µg/mL. The kit offers a sensitivity of 0.129 µg/mL.
Target	LGB
Reactivity	Bovine
UniProt ID	P02754
Application notes	Detection Wavelength: 450 nm
Assay Type	Competitive
Assay Time	1-5h
Range	0.156 µg/mL-10 µg/mL
Sensitivity	0.129 µg/mL
Sample Types	serum, plasma, milk
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)