

## Product Datasheet

### Guinea Pig Alb ELISA Kit (orb339853)

<b>Catalog Number</b>	orb339853
<b>Category</b>	Assays and Kits
<b>Description</b>	The Guinea Pig Alb ELISA Kit detects the target ALB in guinea pig samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 0.156 µg/mL-10 µg/mL. The kit offers a sensitivity of 0.09 µg/mL.
<b>Target</b>	ALB
<b>Reactivity</b>	Guinea pig
<b>UniProt ID</b>	<b>Q6WDN9</b>
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Assay Type</b>	Competitive
<b>Assay Time</b>	1-5h
<b>Range</b>	0.156 µg/mL-10 µg/mL
<b>Sensitivity</b>	0.09 µg/mL
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Note</b>	For research use only
<b>Expiration Date</b>	Please enquire.

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)