

Product Datasheet

Pig Albumin ELISA Kit (orb339884)

Catalog Number	orb339884
Category	Assays and Kits
Description	The Pig Albumin ELISA Kit detects the target ALB in porcine samples. It measures levels in sample types including serum and plasma and urine. The detection range is 0.03 µg/mL-7000 µg/mL. The kit offers a sensitivity of 0.03 µg/mL.
Target	ALB
Reactivity	Porcine
UniProt ID	P08835
Application notes	Detection Wavelength: 450 nm
Assay Type	Competitive
Assay Time	1-5h
Range	0.03 µg/mL-7000 µg/mL
Sensitivity	0.03 µg/mL
Sample Types	serum, plasma, urine
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)