

Product Datasheet

Mouse ACE ELISA Kit (orb340001)

Catalog Number	orb340001
Category	Assays and Kits
Description	The Mouse ACE ELISA Kit detects the target ACE in mouse samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 1.56 ng/mL-100 ng/mL. The kit offers a sensitivity of 0.39 ng/mL.
Target	ACE
Reactivity	Mouse
UniProt ID	P09470
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	1.56 ng/mL-100 ng/mL
Sensitivity	0.39 ng/mL
Sample Types	serum, plasma, tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)