

Product Datasheet

Canine Free Tri-iodothyronine Index ELISA Kit (orb340007)

Catalog Number	orb340007
Category	Assays and Kits
Description	The Canine Free Tri-iodothyronine Index ELISA Kit detects the target FT3 in canine samples. It measures levels in sample types including serum and plasma. The detection range is 2 pmol/L-32 pmol/L. The kit offers a sensitivity of 0.38 pmol/L.
Target	FT3
Reactivity	Canine
Application notes	Detection Wavelength: 450 nm
Assay Type	Competitive
Assay Time	1-5h
Range	2 pmol/L-32 pmol/L
Sensitivity	0.38 pmol/L
Sample Types	serum, plasma
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)