

Product Datasheet

Rat FT4 ELISA Kit (orb340041)

Catalog Number	orb340041
Category	Assays and Kits
Description	The Rat FT4 ELISA Kit detects the target FT4 in rat samples. It measures levels in sample types including serum and plasma. The detection range is 4 pmol/L-64 pmol/L. The kit offers a sensitivity of 1.08 pmol/L.
Target	FT4
Reactivity	Rat
Application notes	Detection Wavelength: 450 nm
Assay Type	Competitive
Assay Time	1-5h
Range	4 pmol/L-64 pmol/L
Sensitivity	1.08 pmol/L
Sample Types	serum, plasma
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)