

Product Datasheet

Bovine Progesterone ELISA Kit (orb340046)

Catalog Number	orb340046
Category	Assays and Kits
Description	The Bovine Progesterone ELISA Kit detects the target PROG in bovine samples. It measures levels in sample types including serum and plasma. The detection range is 1 ng/ml-70 ng/ml. The kit offers a sensitivity of 0.2 ng/mL.
Target	PROG
Reactivity	Bovine
Application notes	Detection Wavelength: 450 nm
Assay Type	Competitive
Assay Time	1-5h
Range	1 ng/ml-70 ng/ml.
Sensitivity	0.2 ng/mL
Sample Types	serum, plasma
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)