

Product Datasheet

Guinea Pig Thyroxine ELISA Kit (orb340075)

Catalog Number	orb340075
Category	Assays and Kits
Description	The Guinea Pig Thyroxine ELISA Kit detects the target T4 in guinea pig samples. It measures levels in sample types including serum and plasma. The detection range is 20 ng/mL-320 ng/mL. The kit offers a sensitivity of 20 ng/mL.
Target	T4
Reactivity	Guinea pig
Application notes	Detection Wavelength: 450 nm
Assay Type	Competitive
Assay Time	1-5h
Range	20 ng/mL-320 ng/mL
Sensitivity	20 ng/mL
Sample Types	serum, plasma
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)