

Product Datasheet

Mouse TNFR1 ELISA Kit (orb340116)

Catalog Number	orb340116
Category	Assays and Kits
Description	The Mouse TNFR1 ELISA Kit detects the target TNFRSF1A in mouse samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 15.625 pg/mL-1000 pg/mL. The kit offers a sensitivity of 1.261 pg/mL.
Target	TNFRSF1A
Reactivity	Mouse
UniProt ID	P25118
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	15.625 pg/mL-1000 pg/mL
Sensitivity	1.261 pg/mL
Sample Types	serum, plasma, tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)